

Optima XE Ultracentrifuge



#

With all the features you expect, along with many advanced functions you won't find anywhere else, the Optima™ XE simplifies use while offering secure loading for efficient, worry-free separations.

Designed as a total system, Optima instruments, rotors and labware are designed, manufactured and tested as an integrated system for optimal performance and safety.

BioSafe	No
Acceleration	10 acceleration profiles
Deceleration	11 deceleration profiles
Timed Run	999 hrs, 59 mins (including hold runs)

Max Capacity Volume Range	- 1500 mL
Temperature Control	± 0.5 °C of set temperature
Temperature Range(s)	Set temperature range: 0 to 40°C in 1°C increments, Ambient temperature range: 10° to 35°C
Display	Full-color LCD touch screen
Item Specifications Referenced	A99833

#